


NOTES: UNLESS OTHERWISE SPECIFIED.



INSTALL CONNECTORS WITH TEFLON TAPE



GASKET ARE TO BE PLACED IN A BAG AND ATTACHED VIA TAPE



REVISIONS

REV	DESCRIPTION	DATE	APP.	REV. BY
0	Init Rel Per ECO 98-0250	7-23-98	KFP	--

ITEM	A	B	C	D	E	F	G	H	J	K	L	M	N	P	R	S	T	U	V	X	Y	Z	AA	AB	AC	AD	AE	AF	AG	PART No.	DESCRIPTION		
7			2				2								2														2	G 284	GASKET, 1/4"		
6			2				2								2														2	N 284	NUT, MALE, 1/4"		
5			2				2								2														2	A 347	1/4 VCR X 1/4T		
4																						2							2	C 137	CONN 1/8T X 1/8NPT BRASS		
4																						2							2	C 66	CONN 1/4T X 1/8NPT BRASS		
4				2				2							2														2	C 64	CONN 1/8T X 1/8NPT SS		
4		2			2				2						2														2	C 966	CONN 6mm X 1/8NPT SS		
4	2		2		2			2		2					2														2	C 67	CONN 1/4T X 1/8NPT SS		
3																2	2	2	2												O 229	O-RING EPDM	
3																2	2	2	2												O 160	O-RING KALREZ	
3								2	2	2	2															2	2	2	2	2	2	O 69	O-RING VITON
3	2	2	2	2																											O 38	O-RING BUNA N STD	
2	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	A-12361	BAFFLE	
1																																B-9729	RECEPTACLE, NYLON
1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	B-13124	RECEPTACLE, SS	

BILL OF MATERIAL

DO NOT SCALE DWG
 TOLERANCE UNLESS OTHERWISE SPECIFIED: ANGULAR ±1/2°
 LINEAR { .X = ±.1
 .XX = ±.02
 .XXX = ±.010

THIS DRAWING IS THE PROPERTY OF TELEDYNE ANALYTICAL INSTRUMENTS AND CONTAINS CONFIDENTIAL INFORMATION. IT IS NOT TO BE COPIED, REPRODUCED OR USED WITHOUT WRITTEN PERMISSION.

TELEDYNE
 ELECTRONIC TECHNOLOGIES
 Analytical Instruments
 An Analyticon Technology Company
 CITY OF INDUSTRY, CALIFORNIA 91748

S/	SIGNATURES	DATE
N/	DRFT: VFig	7/7/98
I/	CHK:	
P/	APPR: K.F. PENSON	7-23-98
O/	ENGR: KARL F. PENSON	
F/	S.O.:	
REFERENCE	CAD I.D. B69583-0	

TITLE	ASSEMBLY, RECEPTACLE (WITH FITTING)
SCALE	NTS
SIM	
SHEET	1 OF 1
MATL.	LISTED
DWG NO.	B-69583
REV	0