





NOTES: UNLESS OTHERWISE SPECIFIED.

1. ALL DIMENSIONS ARE INCHES (±.06) AND MILLIMETERS (±1.) INCHES (mm)

2. LEAVE APPROX. 2" OF CLEARANCE FOR ELECTRICAL CONNECTIONS.

3. PANEL CUT-OUT REQUIRED: 8.75 W. x 11.50 H. (222.) (292.)

REVISIONS				
REV	DESCRIPTION	DATE	APP.	REV. BY
1	INC ECO 97-0386	5/12/97	A.A.	M.V.
2	INC ECO 99-0482	12/9/99	SB	VF



ITEM	QTY	PART NO.	DESCRIPTION
BILL OF MATERIAL			
DO NOT SCALE DWG			
TOLERANCE UNLESS OTHERWISE SPECIFIED: ANGULAR ±1/2°			
LINEAR .XX = ±.1			
.XXX = ±.02			
.LXXX = ±.010			
SIGNATURES		DATE	TITLE
N/	DRFT: JOHN REYES	08/02/96	MDL. 2000 W/ GAS PANEL T/C ANALYZER OUTLINE DIAGRAM
I/	CHK:		
O/	APPR:		
P/ B-67057	ENGR:		
F/	S.O.		
REFERENCE	CAD ID D67056-2		
MATERIAL: NONE			DWG NO. D-67056
SCALE: 1:2			REV: 2
SIM: NONE			SHEET: 1 OF 1