


## Voltage and Frequency Change - Series' 6200, 6400, 9110, 9130 and GFC

### *Part Number Matrix of Voltage / Frequency Dependent Components*

| MODELS   | VOLTAGE              | JP7 (PUMP)        | JP6 (HEATER) | JP2 (AC)  | OPTION    |           |
|--|----------------------|-------------------|--------------|-----------|-----------|-----------|
| <b>6200E</b>   | 100V-120V<br>60HZ    | 042890100         | 040300300    | 040300100 | 050610100 |           |
| <b>6200T</b>   | 100V-120V<br>50HZ    | 042890200         |              |           | 050610200 |           |
| <b>6400E</b>   |                      | 220V-240V<br>50HZ |              | 042890300 | 040300200 | 050610300 |
| <b>6400EH</b>  | 042890400            |                   |              | 050610400 |           |           |
| <b>6400T</b>   | 220V-240V<br>60HZ    |                   |              |           |           |           |
| <b>6400TH</b>  |                      |                   |              |           |           |           |
| MODELS   | VOLTAGE              | PUMP PACK         | JP6 (HEATER) | JP2 (AC)  | OPTION    |           |
| <b>9110E</b><br><b>9110EH</b><br><b>9110T</b><br><b>9110TH</b><br><b>9130E</b><br><b>9130T</b>                       | 100V 60HZ            | 009810600*        | 040300300    | 040300100 | 050610700 |           |
|  | 110V 50HZ            | 009810500*        |              |           |           |           |
|  | 115V 60HZ            | 009810300*        |              |           |           |           |
|  | 230V 50HZ            | 009810400*        |              | 040300200 | 050610900 |           |
|  | 220V-240V<br>50-60HZ | 009810700*        |              |           |           |           |
| MODELS   | VOLTAGE              | PUMP              | HEATER       | JP1 (AC)  | OPTION    |           |
| <b>GFC-7000E</b><br><b>GFC-7000T</b><br><b>GFC-7001E</b><br><b>GFC-7001T</b><br><b>GFC-7002E</b><br><b>GFC-7002T</b> | 100V 50HZ            | 042890200         | N/A          | 038650300 | 043250500 |           |
|  | 100V-120V<br>50HZ    | 042890200         |              | 038650100 | 038650100 | 043250200 |
|  | 100V-120V<br>60HZ    | 042890100         |              |           |           | 043250100 |
|  | 220V-240V<br>50HZ    | 042890400         |              | 038650200 | 038650200 | 043250300 |
|  | 220V-240V<br>60HZ    | 042890300         |              |           |           | 043250400 |

\* - option does NOT include this part

*Relay PCB, All But GFC Series*



*Relay PCB, GFC Series*

